

THE LEADING MANUFACTURER OF PREMIUM
QUALITY CONCENTRIC EMG ELECTRODES



ORYX[®] Healthcare
www.oryxhealthcare.com

The name of our company(Oryx) is inspired by the sharp and unique horns of this rare under extinction animal which indicates the sharpness of our EMG electrodes

"through the animals faced extinction in the early 1970's, when only 6 remained in the wild, the numbers have dramatically increased, thanks to breathing programs and conservation efforts." Ali Faisal Al-Houti










www.nationalgeographic.com/photo-of-the-day/photo/arabian-oryx

CONCENTRIC EMG ELECTRODE

FEATURES

- Precise signal and zero baseline
- Ultra sharp and smooth penetration
- Ergonomic design
- Made of stainless steel with high level of strength



	Plastic Hub Ref number	Metal Hub Ref number	Length (mm)	Diameter (mm)	Surface (mm)	Gauge
	PEN102838	MEN102838	38	0.34	0.02	28
	PEN102828	MEN102828	28	0.34	0.02	28
	PEN102628	MEN102628	28	0.45	0.07	26
	PEN102638	MEN102638	38	0.45	0.07	26
	PEN102650	MEN102650	50	0.45	0.07	26
	PEN102375	MEN102375	75	0.63	0.07	23

Made of stainless steel



easy to open



Ultra sharp and smooth penetration



CONNECTION CABLE



FEATURES

- precise signal transferring
- 5 pin standard connector
(compatible with any EMG machine)
- Highly flexible
- sustainable quality
- reusable
- length : 150 CM



VISION



along side with our current projects, one of our priorities has always been humanitarian activities, for instance, since the very beginning of Russia-Ukraine war we heartfully understood the severance of the situation and despite of all logistical impossibilities, we successfully managed to find solutions and regularly shipped our products (EMG electrodes) to facilitate diagnosing of the injured people

our main vision is to enhance the quality of life by developing innovative projects



Dudullu OSB Mahallesi Des-118 Sk.No:1,Ümraniye
Istanbul, Turkey, 34775

www.oryxhealthcare.com
info@oryxhealthcare.com



[www.linkedin.com
company/oryx-healthcare-ltd](http://www.linkedin.com/company/oryx-healthcare-ltd)



+905312223305